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2. Bild

European Authorized Representative

Certificate of CE registration

QAD 1056

Manufacturer name and address: Hoyotek Biomedical Co. Ltd.
Floor 4, Zone C, Workshop No.1, China civil Aviation science and technology industrialization base No. 225, Jinger Road, Tianjin Airport Economic Zone, China

Product name:	Model:
Corona Virus (COVID-19) Antigen Rapid Test (colloidal gold)	HYT-G01
Influenza A/B/Corona Virus (COVID-19) Antigen Rapid Test (colloidal gold)	HYT-G02
Corona Virus (COVID-19) Combined (IgM/IgG/Neutralizing antibody) Rapid Test (colloidal gold)	HYT-G03

QAdvis EAR as a European Authorized Representative designated by the manufacturer certifies that the products listed above have been notified and filed at the Competent Authority, Swedish Medical Products Agency, as CE-marked In Vitro Diagnostic Medical Devices in accordance with the In Vitro Diagnostic Medical Devices Directive 98/79/EC, article 10.3.

The manufacturer has provided QAdvis EAR with Declaration of Conformity declaring conformance with the requirements in In Vitro Diagnostic Medical Devices Directive 98/79/EC.

Registration data at Swedish Medical Products Agency:

Reference number: 1607413694862
Initial notification to Swedish Medical Products Agency for the products listed above on 2020-12-08.

Date: 2020-12-08

Bing Wu
EAR manager

QAdvis is a member of

QAdvis EAR AB

Address: Ideon Science Park, Scheelevägen 17, SE-223 70 Lund, Sweden
Tel office: +46 8 621 01 05, Email: ear@qadvis.com, Web: www.qadvis.com

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Declaration of conformity

Manufacturer: Hoyotek Biomedical Co., Ltd.
Floor 4, Zone C, Workshop No.1, China civil Aviation science and technology industrialization base No. 225, Jinger Road, Tianjin Airport Economic Zone.

European Representative: QAdvis EAR AB
Ideon Science Park
Scheelevägen 17 SE-223 70 Lund, Sweden

Product Name: Corona Virus (COVID-19) Antigen Rapid Test (Colloidal Gold)
Influenza A/B/Corona Virus (COVID-19) Antigen Rapid Test (Colloidal Gold)
Corona Virus (COVID-19) Combined (IgM/IgG/Neutralizing antibody) Rapid Test (Colloidal Gold)

Product Model: HYT-G01, HYT-G02, HYT-G03

Classification: Other IVD Devices

Conformity assessment Route: IVDD 98/79/EC Annex III

We, Hoyotek Biomedical Co., Ltd hereby declare that the devices mentioned above comply with applicable parts of the Swedish In-Vitro Diagnostic Medical Device Act SFS 1993:584, and the Swedish national legislation LVFS 2001:7, transposing the European In-Vitro Diagnostic Medical Devices Directive, IVDD 98/79/EC.

Verification to: Standard ISO13485:2016, EN ISO14971:2012, EN ISO15223-1:2016, EN ISO18113-1:2011, EN ISO 18113-2:2011, EN ISO18113-3:2011
Related to Directive(s):
98/79/EC (in Vitro Diagnostic Medical Devices)

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【Produktname】

Corona Virus (COVID-19) Antigen Rapid Test (Colloidal Gold)

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- (a) Improper sampling, transport and treatment, and insufficient virus droplets in samples may lead to false negative results.
- (b) Genetic variation in the virus may lead to changes in antigenic determinants, resulting in false negative results, which are more likely to occur with monoclonal antibody reagents.
- (c) For an emergent novel corona virus, the optimal type of sample to be tested and the optimal sampling time after infection (peak viral titer) may not be confirmed. Therefore, multiple sampling at multiple sites in the same patient will reduce the possibility of false negative results.

【Product performance indicator】

- Smooth appearance, solid material attachment, complete contents, complete packaging no damage, clearly identifiable signs, no visible impurities were found in the sample extract.
- The moving speed of sample diluent-10mm per minute.
- Compliance rate of positive quality control products: inspection 5 positive quality control products, P1~P5 Corona Virus (COVID-19) Antigen is required to be positive; the positive internal quality control compliance rate should be 5/5 (+/+).
- Compliance rate of negative quality control products: inspection of 10 negative quality control products, N1~N10 Corona Virus (COVID-19) Antigen is required to be negative; the negative internal quality control compliance rate should be 10/10 (-/-).
- Minimum detection limit: the minimum detection limit of quality control products S1-S5, S1~S4 Corona Virus (COVID-19) Antigen test results should be positive, S5 Corona Virus (COVID-19) Antigen test results should be negative.
- Repeatability: test 2 internal repeatability quality control products, each test for 10 times, test results should be positive.
- The cross reaction of the pathogen: there is no reaction with Flu A virus and Flu B virus.

【Diagnostic performance】

The performance of the kit was tested on samples from patients with suspected COVID-19 taken between March 2020 and January 2021 during daily clinical practice at the Hunan CDC and Shenyang Sixth People's Hospital (China). The samples were taken by qualified personnel. 500 nasopharyngeal swab samples were taken for molecular diagnosis using RT-PCR for the rapid antigen test. The test was shown to have 96 % diagnostic sensitivity and 99 % diagnostic specificity compared with the RT-PCR results.

Clinical Study Results			
Corona Virus (COVID-19) Antigen	PCR Comparator		Sub total
	Positive	Negative	
Positive	96	4	100
Negative	4	396	400
Sub total	100	400	500

Positive Percent Agreement (PPA)= 96/100(96%) (95%CI: 90.07%-98.9%)
 Negative Percent Agreement (NPA)= 396/400 (99%) (95%CI:97.46%-99.73%)

400 nasal swab samples were taken for molecular for molecular diagnosis using RT-PCR for the rapid antigen test, the test was shown to have 95.24 % diagnostic sensitivity and 98.18 % diagnostic specificity compared with the RT-PCR results.

Corona Virus (COVID-19) Antigen	PCR Comparator		Sub total
	Positive	Negative	
Positive	120	5	125
Negative	6	269	275
Sub total	126	274	400

Positive Percent Agreement (PPA)= 120/126(95.24%) (95%CI: 89.92%-98.23%)
 Negative Percent Agreement (NPA)= 269/274 (98.18%) (95%CI:95.79%-99.41%)

1. Interference experiment

The following substances were tested at the concentration shown, and no interference was found.

OTC Throat drop (Halls)	15%	Budesonide	2 mg/mL
OTC Throat drop (Ricola)	15%	Menthol	10 mg/mL
OTC Nasal spray (Afrin)	15%	Mucin	10µg/mL
OTC Nasal spray (VicksSinex)	15%	Mometasone	1 mg/mL

Acetyl salicylic acid	10mg/mL	Ibuprofen	3 mM
Acetaminophen	15mg/mL	Flunisolide	120µg/mL
Afrin Nasal Spray (Oxymetazoline)	4%(v/v)	Tobramycin	80µg/mL
Sodium Cromoglycate	12 mg/mL	Fluticasone	0.4ng/mL
Whole Blood	5%(v/v)	Ritonavir	8.0mg/mL
Chlorpheniramine	5 mg/mL	Abidoc	420mg/mL
maleate Dexamethason	1 mg/mL	Peramivir	1.0mmol/L
Doxycycline hyclate	50µM	Quinine	150µM

2. Cross-reactivity

There was no cross-reaction and interference with the potential cross-reacting microorganisms listed below:

MERS-Coronavirus	Florida/USA-2_Saudi Arabia_2014	5 x 10 ⁴ TCID ₅₀ /ml
Influenza A	H1N1 Denver	2 x 10 ⁷ TCID ₅₀ /ml
	H1N1 WSN/33	1.5 x 10 ⁷ TCID ₅₀ /ml
	H1N1 Pdm-09	2 x 10 ⁷ TCID ₅₀ /ml
	H1N1 New Caledonia	1 x 10 ⁷ TCID ₅₀ /ml
Influenza B	H1N1 New jersey	2 x 10 ⁷ TCID ₅₀ /ml
	Nevada/03/2011	2 x 10 ⁷ TCID ₅₀ /ml
	B/Lea/40	5 x 10 ⁷ TCID ₅₀ /ml
Human Coronavirus	B/Taiwan/2/62	1 x 10 ⁷ TCID ₅₀ /ml
	229E	1 x 10 ⁷ TCID ₅₀ /ml
	OC43	1 x 10 ⁷ TCID ₅₀ /ml
Respiratory syncytial virus	NL63	1 x 10 ⁷ TCID ₅₀ /ml
	Type A	1 x 10 ⁷ TCID ₅₀ /ml
	Type B	1 x 10 ⁷ TCID ₅₀ /ml
Human Metapneumovirus (hMPV)	hMPV 3 Type B1 / Peru2-2002	1 x 10 ⁷ TCID ₅₀ /ml
	hMPV 16 Type A1 / IA10-2003	1 x 10 ⁷ TCID ₅₀ /ml
Parainfluenza virus	Type 1	1 x 10 ⁷ TCID ₅₀ /ml
	Type 2	1 x 10 ⁷ TCID ₅₀ /ml
	Type 3	1 x 10 ⁷ TCID ₅₀ /ml
	Type 4A	1 x 10 ⁷ TCID ₅₀ /ml
Rhinovirus	A16	1 x 10 ⁷ TCID ₅₀ /ml
	Type B42	1 x 10 ⁷ TCID ₅₀ /ml
Enterovirus	Type 68	1 x 10 ⁷ TCID ₅₀ /ml
	(09/2014 isolate 4)	1 x 10 ⁷ TCID ₅₀ /ml
Mycobacterium tuberculosis	K	1 x 10 ⁷ TCID ₅₀ /ml
	Erdman	1 x 10 ⁷ TCID ₅₀ /ml
	HN878	1 x 10 ⁷ TCID ₅₀ /ml
	CDC1551	1 x 10 ⁷ TCID ₅₀ /ml
	H37Rv	1 x 10 ⁷ TCID ₅₀ /ml

Novel Coronavirus Monoclonal Antibody is a protein that recognizes nucleocapsid protein and can detect genetic variants of strains.

【Cautions】

- The test should not be done if the condition of kit and samples are not restored to 18-30°C, as it will affect the accuracy of the results.
- The positive samples obtained by the rapid test should be confirmed by other methods.
- The rapid test should be sealed and kept in dry place. The test cassette should be tested as soon as possible after removal from the packaging. Avoid placing it in the air for too long due to moisture.
- The deepness of the test line color is not necessarily associated with the titer of the antigen in the sample. The interpretation of the results after 15 minutes are invalid.
- When COVID-19 antigen content in the sample is very high, the C-line zone may be weakened, which is a normal phenomenon.
- The results of rapid test are only for clinical reference and should not be the only basis for clinical diagnosis and treatment.
- Waste samples and test should be treated as potential infectious agents.
- The appearing time of the Control Line should not be taken as the time basis for judging the results of test line. The color rendering results should be observed and judged within a time limit of 13-15 minutes.

- The rapid test is only used for in vitro diagnosis.
- This product must be operated by professionally trained personnel, such as medical staff with clinical experience.

Hoyotek Biomedical Co., Ltd
 Floor 4, Zone C, Workshop No.1, China civil Aviation science and technology industrialization base No. 225, Jinger Road, Tianjin Airport Economic Zone. 300308 Tianjin China.

QAdvis EAR AB
 Ideon Science Park
 Scheelevägen 17, SE-223 70 Lund, Sweden

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	Do not reuse		For in vitro diagnostic use only
	Store between 2-30°C		Consult instructions for use
	Use by		Lot number
	Do not use if package is damaged		Contains sufficient for <n> tests
	Keep away from sunlight		Keep dry
	Manufacturing date		Manufacturer
	Authorized representative in the European Community		

Version Number: 02
 Effective Date: June 3rd, 2021

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Hoyotek Biomedical Co., Ltd.	Doc No.	HYT-G01-009
	Rev.:	00
	Page:	1-11

Clinical evaluation report

Company Name: Hoyotek Biomedical Co., Ltd.

Packing specification: 30 tests/box

Company Address: Floor 4, Zone C, Workshop No.1, China civil Aviation science and technology
industrialization base No. 225, Jinger Road, Tianjin Airport Economic Zone.

Product Name: Corona Virus (COVID-19) antigen Rapid Test (Colloidal gold)

Clinical evaluation place: the National CDC of Hunan and Shenyang Sixth People's Hospital

Start date: Mar 10, 2020

End date: January 19, 2021

Hoyotek Biomedical Co., Ltd.

Hoyotek Biomedical Co., Ltd.	Doc No.	HYT-G01-009
	Rev.:	00
	Page:	2-11

Introduction

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped non-segmented positive-sense RNA virus. It is the cause of coronavirus disease (COVID-19), which is contagious in humans. SARS-CoV-2 has several structural proteins including spike (S), envelope (E), membrane (M) and nucleocapsid (N). The antigen is generally detectable in upper respiratory samples during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease Perform clinical evaluation of the Corona Virus (COVID-19) antigen Rapid Test (Colloidal gold) to determine its diagnostic sensitivity and specificity.

Detection principle

The product is based on the principle of sandwich and colloidal gold immunochromatography, the nitrocellulose membrane Test Zone is pre-coated with mouse anti-Corona Virus (COVID-19) monoclonal antibody, the Control Zone is pre-coated with goat anti-mouse polyclonal antibody, the gold conjugation pad is pre-coated with colloidal gold labeled mouse anti-Corona Virus (COVID-19) monoclonal antibody. When testing positive samples, COVID-19 Antigen the sample will be combined with colloidal gold (Au) labeled mouse anti-Corona Virus (COVID-19) monoclonal antibody and form into immune complexes (Au-mouse anti-Corona Virus (COVID-19) monoclonal antibody-[COVID-19-(Antigen)]), the complexes will move forward inside the nitrocellulose membrane by chromatography effect. When reaching Test Zone, the complexes will be combined with mouse anti-Corona Virus (COVID-19) monoclonal antibody, and form into “(Au-mouse anti-Corona Virus (COVID-19) monoclonal antibody-[COVID-19-(Antigen)]-[mouse anti-Corona Virus (COVID-19) monoclonal antibody]”, thus agglutination color appears; The residual colloidal gold labeled mouse anti-Corona Virus (COVID-19) monoclonal antibody will be combined with goat anti-mouse polyclonal antibody at the Control Zone, and produce color under agglutination. When testing negative samples, there’s no COVID-19 Antigen in the samples, no immune complexes will be formed and color only appears in the Control line.

Purpose

Perform clinical evaluation of the Corona Virus (COVID-19) antigen Rapid Test (Colloidal gold) to determine its diagnostic sensitivity and specificity.

Testing management

During the trial, the main investigator is responsible for the coordination and management of the entire clinical trial, and the main participants are responsible for the main trial work. During the clinical trial, the main researcher supervises the quality control of the testing laboratory. Any problems found in the test must be contacted with the main researcher in time and appropriate measures should be taken. The final test results are statistically analyzed by the person in charge of statistics, and the main investigator confirmed and wrote the report.

Hoyotek Biomedical Co., Ltd.

Hoyotek Biomedical Co., Ltd.	Doc No.	HYT-G01-009
	Rev.:	00
	Page:	3-11

Methods

Synchronous blind test and methodological comparison design. The Nasopharyngeal swab and oropharyngeal swab samples and saliva samples were collected by hospital professional medical staff in accordance with the sampling methods of Corona Virus (COVID-19) antigen Rapid Test (Colloidal gold) and Sansure Biotech Inc novel coronavirus detection kit. The Nasopharyngeal swab and oropharyngeal swab samples are blindly numbered and grouped by the National CDC of Hunan and Shenyang Sixth People's Hospital. Nasopharyngeal swab samples are divided into one group, saliva samples are divided into one group, oropharyngeal swab samples are divided into one group, and then tested by laboratory inspectors

Discussion and Conclusion

Results

In this clinical trial, Nasopharyngeal swab specimens were obtained from e National CDC of Hunan and Shenyang Sixth People's Hospital and tested with the SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold) and the comparator device Sansure Biotech Inc novel coronavirus detection kit produced. Statistical analysis was performed to calculate the positive agreement rate and negative agreement rate.

In this study, The Corona Virus (COVID-19) antigen Rapid Test (Colloidal gold) produced by Hoyotek has a diagnostic sensitivity of 96% for detecting Nasopharyngeal swabs, 93% for detecting oropharyngeal swab, 91% for detecting Saliva. The diagnostic specificity of 99% for detecting Nasopharyngeal swabs, 97.5% for detecting oropharyngeal swab, and 97.3% for detecting Saliva.

The results showed that the investigational device, the Antigen Rapid Test Kit (Colloidal Gold), meets the needs of clinical testing.

Main Content

General design

This test uses a synchronous blind test and methodological comparison design. In order to eliminate the possible impact of the subjective biases and personal preferences of researchers on the test results during the clinical trial process, this test uses a blind test. That is, the test personnel in this test do not know the specific information of the sample, and the clinical information of the sample may not be released until the end of the test. After the samples were enrolled, the samples were coded by the blind editor authorized by the clinical trial, in which the blind editor was not involved in the test operation of the clinical trial. Testing personnel shall test the coded sample according to the reagent test specification. In the process of test operation, clinical test researchers should strictly follow the requirements of the product specification for test operation and interpretation check, and the results obtained in the test process should be truthfully recorded in the data collection table.

For the detection of SARS-CoV-2 Antigen, the Nasopharyngeal swab and r oropharyngeal swab samples were collected by hospital professional medical staff in accordance with the sampling methods of Corona

Hoyotek Biomedical Co., Ltd.

Hoyotek Biomedical Co., Ltd.	Doc No.	HYT-G01-009
	Rev.:	00
	Page:	4-11

Virus (COVID-19) antigen Rapid Test (Colloidal gold) and Sansure Biotech Inc novel coronavirus detection kit. The nasopharyngeal swab and oropharyngeal swab samples are blindly numbered and grouped by the National CDC of Hunan and Shenyang Sixth People's Hospital. Nasopharyngeal swab samples are divided into one group, saliva samples are divided into one group, oropharyngeal swab samples are divided into one group, and then tested by laboratory inspectors.

Measures to reduce and avoid bias

Subjects were screened strictly according to the blind grouping of the clinical trial protocol to reduce the selection bias.

Prior to the start of the trial, the lab operators were trained to correctly perform the tests and follow the trial protocol.

Clinical sample related requirements

1) DOs and DON'Ts of Sample Collection

- Do test sample immediately
- Use only swabs provided with the kit.

2) Sample storage

- Specimen Transport and Storage
- Freshly collected specimens should be processed within 1 hour, or stored at the manufacturer's recommended condition before testing.
- It is essential that correct specimen collection and preparation methods be followed.

Clinical sample selection

1) Inclusion criteria

Sample inclusion criteria: the sample should be a sample with clearly recorded source, including different age, gender and other factors. The collection and treatment of samples are in accordance with the reagent specification or relevant regulations. Sample information should be complete, including age, sex, sample collection date, clinical diagnosis such as confirmation or exclusion of SARS-CoV-2 infection.

2) Exclusion criteria

- Samples that are unable to complete the test process human factors (sample contamination during operation).
- Samples were contaminated with bacteria or/and nosebleed.
- Samples went through too many freeze-thaw cycles.
- Samples not kept at the requirement conditions.

Hoyotek Biomedical Co., Ltd.	Doc No.	HYT-G01-009
	Rev.:	00
	Page:	5-11

➤ Samples that are unable to yield control signals with PCR reagent.

Quality control

Definition

Quality control is defined as the operation of techniques and activities, such as monitoring, under the quality assurance system to verify that the research quality meets the requirements. Quality control must be applied at every stage of data processing to ensure that all data is trusted and properly located.

1) Study monitoring

During the outbreak, authorized and qualified inspectors will conduct regular remote 6 primary data checks according to the monitoring plan to verify compliance with protocols and regulations and assist investigators.

2) Laboratory quality control

The laboratory of the testing shall establish a unified test index, standard operating procedures and quality control procedures.

3) Quality control of reagent testing process In each test, the control line shall have red strip (qualified quality control). If the control line does not have red strip (unqualified quality control), the cause shall be found out and retested until the quality control result is qualified, so as to ensure the reliability and stability of the system.

4) Qualification of researchers The researchers participating in the clinical trial must have the specialty, qualification and ability of the clinical trial, and pass the qualification examination. The personnel requirements should be relatively fixed.

Reagents and instruments for clinical research

The information of reagents for test is shown in Table 1:

Table 1 Reagent Information

S/N	Sample Name	Size	Lot No.	Company
1	Corona Virus (COVID-19) Antigen Rapid Test (Colloidal Gold)	30/box	20201001	Hoyotek

Determination of results: Operate according to the instructions. After the test, observe the color reaction of colloidal gold on the T line and C line to determine the test result of the new coronavirus (COVID-19) antigen.

The nucleic acid background of all swab samples are the results of the new coronavirus nucleic acid detection reagent produced by Sansure Biotech Inc.

Clinical Trial Results and Analysis

Hoyotek Biomedical Co., Ltd.

Hoyotek Biomedical Co., Ltd.	Doc No.	HYT-G01-009
	Rev.:	00
	Page:	6-11

Clinical diagnostic sensitivity of reagents

Table 2 Results of 100 cases of samples tested positive for COVID-19

	Nasopharyngeal swabs Number of positive	oropharyngeal swab Number of positive	Saliva Number of positive
Hoyotek antigen detection results	96	93	91

Test result: The Corona Virus (COVID-19) antigen Rapid Test (Colloidal gold) produced by Hoyotek have a diagnostic sensitivity of 96% for detecting Nasopharyngeal swabs, 93% for detecting oropharyngeal swab, and 91% for detecting Saliva.

Diagnostic Specificity

Table 6 Results of 400 Nasopharyngeal /oropharyngeal swab samples with negative COVID-19 nucleic acid test

	Nasopharyngeal swabs Number of positive	oropharyngeal swab Number of positive
Hoyotek antigen detection results	396	390

The Corona Virus (COVID-19) antigen Rapid Test (Colloidal gold) produced by Hoyotek detects 400 cases nucleic acid negative Nasopharyngeal swabs with a diagnostic specificity of 99%, oropharyngeal swabs with a diagnostic specificity of 97.5%

Table 7 Results of 151 saliva samples with negative COVID-19 nucleic acid test

	Number of negative	Number of positive
Hoyotek antigen detection results	147	4

The Corona Virus (COVID-19) antigen Rapid Test (Colloidal gold) produced by Hoyotek detects 151 cases nucleic acid negative saliva with a diagnostic specificity of 97.3%.

Conclusion:

This clinical trial has performed a full analysis of the experimental reagents through methodological comparisons, and the results all meet the criteria for clinical evaluation. All the results showed that SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold) meet the needs of clinical test.

1) Appendix 1

S/N	age	Gender	Types	collecti on date	ORF-Lab Scope- Fam CT Value	N gene-ROX CT Value	internal standard - HEX CT	clinical symptoms
-----	-----	--------	-------	---------------------	--------------------------------	------------------------	----------------------------------	----------------------

Hoyotek Biomedical Co., Ltd.

							Value	
1	63	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	22.41	30.02	27.74	Serious
2	29	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	33.39	31.85	26.93	slight
3	25	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	31.24	31.72	29.08	slight
4	62	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	37.85	31.51	29.75	slight
5	77	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	31.81	33.13	25.07	slight
6	36	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	40.26	33.05	28.06	slight
7	55	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	20.79	35.67	26.67	slight
8	73	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	36.33	30.90	26.88	Serious
9	25	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	39.29	35.33	28.45	slight
10	20	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	33.26	31.35	29.31	slight
11	68	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	23.95	37.44	27.03	medium
12	69	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/11	23.72	31.68	28.77	slight
13	25	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/12	31.27	37.58	26.62	slight
14	50	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/12	28.01	33.17	28.84	slight
15	39	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/12	23.73	33.01	27.06	slight
16	61	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/12	39.00	31.71	29.08	Serious
17	28	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/13	32.76	38.51	25.22	medium
18	21	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/13	35.23	35.20	26.15	medium
19	43	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/13	22.53	37.71	29.12	slight
20	36	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/13	32.78	33.34	25.99	slight
21	52	Man	Nasopharyngeal	2020/3/	20.99	32.77	29.93	slight

			swabs/Oropharyngeal /saliva	13				
22	32	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 13	26.50	38.17	26.06	slight
23	25	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 13	31.75	34.33	27.51	medium
24	74	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 13	26.03	37.82	26.08	Serious
25	38	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 13	36.66	36.54	26.55	slight
26	21	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 13	31.31	31.55	25.17	slight
27	26	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 13	22.33	31.93	27.44	slight
28	26	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 13	36.41	37.92	27.28	slight
29	77	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	22.10	32.95	27.60	slight
30	29	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	21.21	30.47	28.35	slight
31	71	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	28.19	34.37	29.11	Serious
32	46	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	30.30	30.73	28.10	slight
33	57	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	22.76	31.04	25.90	slight
34	47	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	25.78	31.86	27.84	slight
35	52	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	30.04	34.92	25.88	slight
36	64	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	34.65	36.35	27.29	slight
37	48	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	39.48	38.91	27.63	slight
38	70	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	30.05	33.50	27.64	slight
39	61	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	27.18	33.45	28.48	medium
40	37	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 14	34.31	35.49	28.83	medium
41	25	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	23.95	38.06	29.85	medium
42	65	Woman	Nasopharyngeal	2020/3/	21.48	32.03	29.12	slight

Hoyotek Biomedical Co., Ltd.	Doc No.	HYT-G01-009
	Rev.:	00
	Page:	9-11

			swabs/Oropharyngeal /saliva	15				
43	28	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	36.04	37.41	27.41	medium
44	56	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	37.45	36.29	29.73	slight
45	37	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	26.58	30.73	27.56	slight
46	60	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	28.40	31.13	25.30	medium
47	59	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	33.43	31.60	25.87	slight
48	34	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	28.21	31.55	26.89	slight
49	52	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	31.12	30.08	26.24	slight
50	67	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 15	27.30	32.33	27.28	slight
51	40	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 16	36.97	36.03	29.81	slight
52	40	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 16	26.68	38.15	27.21	slight
53	32	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 16	23.84	32.68	25.27	slight
54	31	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 18	23.57	31.87	29.80	slight
55	22	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 18	33.51	35.89	26.22	slight
56	57	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 18	38.63	30.13	29.33	slight
57	64	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 18	29.22	36.18	29.48	slight
58	46	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 18	34.23	31.87	26.28	medium
59	72	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 18	36.49	33.51	27.48	slight
60	27	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 18	38.41	32.73	26.08	medium
61	43	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 19	28.00	34.32	26.70	slight
62	35	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 19	26.50	35.62	27.07	slight
63	56	Woman	Nasopharyngeal	2020/3/	33.93	30.81	28.58	slight

Hoyotek Biomedical Co., Ltd.

			swabs/Oropharyngeal /saliva	19				
64	28	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 19	35.01	36.48	29.61	slight
65	23	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 20	20.45	36.22	25.21	medium
66	48	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 20	35.26	37.81	29.78	slight
67	31	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 20	21.84	36.38	25.69	slight
68	21	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 20	22.96	30.29	26.02	slight
69	28	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 20	26.19	38.33	25.46	slight
70	71	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 20	26.16	38.32	29.06	slight
71	46	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 20	32.88	37.66	28.91	slight
72	25	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 20	31.88	38.06	28.06	slight
73	51	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 21	23.99	36.61	27.81	medium
74	31	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	24.22	31.53	28.28	slight
75	37	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	22.18	33.66	25.19	slight
76	29	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	36.01	31.97	29.95	slight
77	20	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	39.52	30.48	25.91	slight
78	75	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	21.72	31.29	25.23	slight
79	62	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	36.51	34.46	28.05	slight
80	35	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	37.66	38.92	26.56	medium
81	60	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	33.17	35.72	28.40	slight
82	59	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	20.30	37.31	27.77	slight
83	22	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 22	39.77	38.09	29.43	slight
84	26	Woman	Nasopharyngeal	2020/3/	37.71	32.68	26.56	medium

Hoyotek Biomedical Co., Ltd.	Doc No.	HYT-G01-009
	Rev.:	00
	Page:	11-11

			swabs/Oropharyngeal /saliva	23				
85	37	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	33.48	34.11	25.15	slight
86	35	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	39.43	37.81	27.50	slight
87	48	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	21.70	30.25	29.72	slight
88	27	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	27.75	31.96	25.28	slight
89	76	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	33.36	34.45	29.02	Serious
90	71	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	20.49	31.05	29.85	slight
91	35	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	20.06	30.28	27.85	slight
92	50	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	36.14	31.59	27.12	slight
93	24	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	34.77	30.46	28.85	slight
94	47	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	34.84	36.37	26.13	medium
95	60	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	20.69	35.48	28.94	medium
96	46	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	40.00	35.11	25.26	slight
97	33	Man	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	29.97	31.15	29.93	slight
98	31	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	37.78	35.38	27.53	medium
99	53	Woman	Nasopharyngeal swabs/Oropharyngeal /saliva	2020/3/ 23	29.20	38.22	26.89	slight
100	30	Woman	Oropharyngeal swabs	2020/3/ 23	26.96	35.66	27.90	slight